

Transformation of *P. pastoris*

Introduction

In Rot sind andere Werte.

Materials

›

- › Electrocompetent *P. pastoris* cells
- › 7-10µg of linearized pPICZα DNA
- › 0.2 cm electroporation cuvette
- › Sterile ice-cold (0°C) 1M sorbitol
- › Electroporation device
- › YPDS + Zeocin plates

›

Procedure

1. Mix 80 µl of the cells from Step 6 (Kompetenz *P. pastoris*) with 5–10 µg of linearized pPICZα DNA (in 5–10 µl sterile water) and transfer them to an ice-cold (0°C) 0.2 cm electroporation cuvette.
2. Incubate the cuvette with the cells on ice for 5 minutes.
3. Pulse the cells according to the parameters for yeast (*Saccharomyces cerevisiae*) as suggested by the manufacturer of the specific electroporation device being used. (or 1.5 kV, 200 Ω, 25 µF)
4. Immediately add 1 ml of ice-cold 1 M sorbitol to the cuvette. Transfer the cuvette contents to a sterile 15 ml tube.
5. Let the tube incubate at 30°C without shaking for 1 to 2 hours. (longer incubation will increase transformation efficiency)
6. Spread 50–200 µl each on separate, labeled YPDS plates containing the appropriate concentration of corresponding antibiotic. (containing 100 µg/ml of Zeocin)
7. Incubate plates for 2 to 3 days (4-5 days) at 30°C until colonies form.
8. Pick 10–20 colonies and purify (streak for single colonies) on fresh YPD or YPDS plates containing the appropriate concentration of corresponding antibiotic.